ENG

First things first

Before using your menstrual cup, check that the 4 holes under the rim of the cup are well pierced and ensure you sterilise it properly, to do so:

- 1) Wash your hands.
- 2) Wash the cup with warm water and mild soap. Rinse thoroughly after.
- 3) Immerse the cup fully in a pan of water and boil it for 5 to 7 minutes.
- 4) Allow to cool and remove it.

Your cup is now sterilised and ready to be used!

Finding your fold

Everyone's unique and when it comes to folding menstrual cups ready for insertion, there are several different methods and we've explained a few below with the first diagrams at the bottom of the leaflet. Remember, it's totally normal if it takes a few tries before you find the fold that works for you.

The C-Fold

The clue's in the title: this one makes your cup look like the letter C.

- 1) Pinch your cup in the middle so it flattens.
- 2) Pinch the outside edges together to make the C-shape.

The Punch-Down Fold

This is a little more complicated than the C-fold but is a very popular method.

- 1) Holding the cup in one hand, use your index finger to push the rim of the cup to the very base.
- 2) Keeping that finger on the base, pinch the sides together to make a narrow opening at the top of the cup.

The Triangle/7 Fold

Okay, so this one sounds more complicated than it is – you do not have to fold the cup 7 times!

- 1) Using both hands, pinch and flatten the cup.
- 2) Take one corner and fold it to the base of the cup, making a triangle.

OR

- 3) Take that one corner and fold it halfway to the base, making a '7' shape'.

Top tip: Folding all the way to make the triangle shape creates a smaller tip/entry point, which is ideal for beginners.

Inserting your cup

Now you've tried the different folds, it's time to insert your folded menstrual cup. You'll need to wash your hands, take a comfortable position (for example, sitting on the toilet seat or squatting) and take a few deep breaths to relax. The more relaxed you are, the easier it is to insert! If you're trying cups for the first time, try using it on your heaviest flow days (usually day 2 or 3 of your period) as this will help with insertion.

Top tip: You may find wetting your cup with a little water helps lubricate it, making insertion that little bit easier. Don't use any oil-based or petroleum jelly lubricants on your cup.

Before we begin, it's important to note a menstrual cup is not designed to sit as high in the vagina as a tampon. In fact, it is designed to sit in the lower part of the vagina – you can see the difference in diagram 2 (Cup Position) below. If it's in the right place it shouldn't be uncomfortable, and you shouldn't be able to feel it inside.

Now you have the cup folded in your preferred method, you must hold onto it firmly, as close as possible to the base of the cup, then using diagram 3 (Inserting your Cup) below:

- 1) Position the top of the cup at the opening of the vagina.
- 2) Aiming towards your lower back, not straight up.
- 3) Insert slowly.

Once you've got the cup in the right place, you may find it hasn't fully unfolded open. It's important the cup 'pops' open when inside as this creates a seal against the walls of your vagina to prevent leaks. You can open the cup fully by pinching the base tightly and twisting the cup around 360°.

Top tip: You can run your finger around the base of the cup to check it's fully open. If you feel an unopened section, gently press right next to it until it pops open.

1. Finding your fold / Trouvez le pli qui vous convient**The C-Fold / Le pli en C****Removing your menstrual cup**

Your cup should be emptied, cleaned and reinserted every 4 to 8 hours. Before removing, wash your hands, get in a comfortable position, and take a few deep breaths.

Before removing the cup completely, you must release the suction that the opened cup has created, to do so:

- Pinch the base of the cup to release the suction and 'wiggle' the cup down to remove.
- When you get near the opening of the vagina, carefully tilt the cup and remove one side of the rim first and then the other. This should help you contain the menstrual fluid within the cup.
- If you cannot reach the base of your cup to remove it, give the stem a gentle pull whilst pushing down with your vaginal muscles until you can pinch the base of the cup.
- Once fully removed:
 - Empty the cup's contents into the toilet and flush.
 - Clean your cup using the instructions below.

Cleaning your cup

When you remove your cup mid-period, wash the cup with cold water and mild soap. Ensure the suction holes are clear and then reinsert.

Top tip: If you're out and about, take bottled water with you and give the cup a thorough rinse with this. Next time you remove your cup, make sure you give it a clean with mild soap and water.

At the end of your period, you'll need to resterilise your cup (as you did before your first use) and dry it thoroughly. Then keep it stored in your cotton bag, in a cool, dry place until your next period.

Word to the wise

It's important to know that menstrual cups can be associated with Toxic Shock Syndrome (TSS). TSS is a rare but serious illness that may be fatal. To reduce the risk of TSS we recommend emptying and cleaning your cup at least every 4 to 8 hours and to resterilise it between periods. Symptoms of TSS can develop quickly (and may seem a bit like flu to begin with), they include: a sudden high fever (usually over 39°C), vomiting, diarrhoea, sunburn-like rash, sore throat, muscle aches, dizziness or fainting. If you experience any of these symptoms during your period, remove your cup and seek medical attention. Ensure you tell whoever treats you that you have been using internal period products. Do not use a menstrual cup if you have ever had TSS.

Three golden rules:

- Don't share your menstrual cup with anyone else.
- A menstrual cup does not prevent pregnancy nor sexually transmitted diseases.
- Any pain, discomfort or difficulty urinating, please seek medical advice.

Also, you can use a cup if you have a coil/IUD fitted. However, we recommend that you speak to your doctor or gynaecologist before usage as, in rare cases, cups can dislodge IUDs.

Customise your cup

After a few uses of your cup, you may think the stem needs adjusting for a more comfortable fit. Our cup has a particularly long stem, which is good for those with a high cervix. However, if you find the cup is in a comfortable position but the stem sits out of your vagina and feels uncomfortable, it may need a trim.


Do not trim the stem while the cup is inside the vagina. Use clean scissors to trim and take a little off at a time, trying it for comfort and trimming more if needed.

Top tip: When you become a menstrual cup pro, you may find you want to trim the stem off completely, but we certainly wouldn't be this bold if it's your first time!

Lifetime of your cup

If cared for as per this leaflet, your menstrual cup should last you many years but check it regularly for damage and wear. A little bit of staining is normal but if your cup splits or becomes sticky then it needs replacing. Your old cup should be disposed of with general household waste.

To find out more about the Lil-Lets products that are available to you, please visit www.lil-lets.com.

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lil-lets®
reusable
menstrual cup

Welcome to your Lil-Lets menstrual cup – great choice! Our cup is made of medical grade silicone and is designed to be super soft and flexible, providing ultimate comfort and making it ideal for beginners.

Menstrual cups are designed to collect menstrual fluid rather than absorb it. You can wear your cup for 4-8 hours before emptying, cleaning, and reinserting.

Lil-Lets reusable menstrual cup
has been specially designed to be ideal for beginners:

Discover other products in the Lil-Lets reusable and organic ranges at www.lil-lets.com

Made in France and designed by Si-Bell for:

Lil-Lets UK Ltd,
PO BOX 14568, Solihull,
B91 9LN, United Kingdom;

Lil-Lets EU,

Unit 3D North Point House, New Mallow Road,
Cork, T23 AT2P, Ireland;

Premier FMCG (Pty) Ltd,
Building 5, Maxwell Office Park, Magwa Crescent West, Waterfall City, 2090, South Africa.


ZULU Siyakwamukela enkomishini yakho yokuya esikhathini yakwaLil-Lets – ukhethe kahle! Inkomishi yethu yenzinwe ngesilikhoni yohlolo lwezokwelapha futhi yenzelwe ukuthi ithambhe okuhulu futhi iwayi kuzukupotopota, ukunikeze ukunetheteka lwaselikhathi zoku esikhathini zenvelo ukucouqa ukucheteki. Izwinkomishi zoku esikhathini kuzukupola zilumunha. Ungayapha inkomishi yakho amahora ama-4 kuya kwayi-8 ngaphambi kokyuthulula, ukuyihlanza nokuphinda uyifaka kabusha.

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